
BiPAP[®] Disposable Circuit Reorder #312107 Assembly Instructions

The BiPAP Disposable Circuit is made up of Circuit Tubing and a Proximal Pressure Line Assembly.

Caution: US federal law restricts this device to sale by or on the order of a physician.

Warnings and Cautions

- The exhaust vent on the Exhalation Port is designed to exhaust CO₂ from the patient circuit. Continuous flow is required for safe operation. **Do not block or otherwise try to seal the exhaust vent on the exhalation port.**
- Before instituting mechanical ventilation, confirm the proper connection of the circuit, and verify ventilator operation. Patient, ventilator and circuit must be monitored on a regular basis per established hospital standards of care.
- Do not pull or stretch the circuit tubing. This could result in circuit leakage.
- This circuit is not recommended for use with a heated humidifier since it does not allow for the use of water traps.

Assembly Instructions

1. CONNECT THE CIRCUIT TUBING.

Connect a low-resistance, mainflow bacteria filter on the outlet port of the BiPAP unit. Connect the cuffed end of the 6' Tubing to the outlet port of the filter.


2. CONNECT THE PROXIMAL PRESSURE LINE ASSEMBLY. (Refer to illustrations below.)

- BiPAP S-D or S/T-D Units with Airway Pressure Monitor (APM) (Figure 1)
 - Connect the 24" length of line to the Pressure Port on the front of the BiPAP unit.
 - Connect the 2" length of line to the Airway Connection Port on the rear of the APM.
- BiPAP S-D or S/T-D Units without Airway Pressure Monitor (APM) (Figure 2)
 - When the Airway Pressure Monitor is not used, only the 84" length of proximal pressure line tubing is needed. Disconnect the 2" and 24" length from the assembly by cutting the 84" tubing near the Tee Connector. Discard the Tee Connector, 2" length and 24" length of tubing.
 - Connect the 84" length to the Pressure Port on the front of the BiPAP unit.
- BiPAP S or S/T Units with Airway Pressure Monitor (APM) (Figure 3)
 - When a non-diagnostic unit is used (BiPAP S unit or BiPAP S/T unit), only the 84" length of proximal pressure line tubing is needed. Disconnect the 2" and 24" length from the assembly by cutting the 84" tubing near the Tee Connector. Discard the Tee Connector, 2" length and 24" length of tubing.
 - Connect the 84" length of line to the Airway Connection Port on the rear of the APM.

3. CONNECT THE EXHALATION PORT END OF THE TUBING TO THE MASK'S ROTATING CONNECTOR.

The vent on the Exhalation Port should be positioned away from the patient's face.

THIS CIRCUIT IS DESIGNED FOR SINGLE PATIENT USE. IT IS NOT INTENDED TO BE CLEANED OR STERILIZED. The length of time a circuit is used should be determined by the appropriate personnel, based on accepted infection control procedures.

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